

# Serum dextrose agar (SDA)

## Détails

For 1 liter of SDA:

|  |    |  |
|--|----|--|
| 33 g of Difco Tryptose blood agar base | or | 10 g of Tryptose<br>3 g of Beef extract<br>5 g of NaCl<br>15 g of Agar |
|--|----|--|

Complete to 950 ml with water.  
Autoclave.

Add 50 ml of sterile horse serum.  
Add 50 ml g of sterile glucose 20%.

## Reference:

Corbel MJ, Thomas EL. 1980.

The Brucella-phages: their properties, characterization and applications. Booklet 2266, Ministry of Agriculture, Fisheries and Food, Pinner, Middlesex, UK, 106 p

Last modifications: 2011-02-24